

## X-ray Photoelectron Spectrometer

### Introduction

**X-ray Photoelectron Spectroscopy** or **Electron Spectroscopy for Chemical Analysis (ESCA)** is a powerful technique to determine the quantitative elemental composition of sample surfaces. XPS can be used to analyze inorganic compounds, metal alloys, semiconductors, polymers, glasses, ceramics, teeth, bones, medical implants, bio-materials, industrial samples, and many others.

### Specifications

<b>Analyzer</b>	<b>180° double-focusing hemispherical analyzer with 128 channel detectors</b>
<b>X-ray source</b>	Al-K $\alpha$ (1486.8 eV) micro-focused monochromator with variable spot size (50-400 $\mu$ m)
<b>Ion Gun</b>	Ar+ Energy range 100-4000 eV
<b>Charge Compensation</b>	Dual-beam source
<b>Sample Stage</b>	60 x 60 mm <sup>2</sup> ; 5 axis auto-z can accommodate samples of maximum size 60 x 60 mm <sup>2</sup> with 15 mm thickness. However, ~5 mm <sup>2</sup> area with < 1 mm thickness is preferred for solid samples.



## Samples

- Solids, crystals, thin films, epitaxial layers, polymers, powders - tightly packed pellets of 6 mm to 10 mm diameter, glasses, ceramics, teeth, bones, medical implants, and biomaterials.
- **Powders / Nanomaterials pressed in Indium metal foil/ the ball is recommended.** If sample preparation is needed user may need to pay additionally for sample preparation.
- Samples should be clean and dry, the samples should not be outgassing in nature, and should not degrade under high vacuum
- Samples composed of Na, K, S, P, F, Zn, Pb, Bi, I, Se, As, Te, and Hg are not suitable for Depth profile analysis
- No liquid and moisture samples shall be done
- If any samples contain **Iodine** and **Sulphur**, XPS measurements cannot be done. However, Sulphur bonded strongly i.e. GaS, MoS<sub>2</sub>, etc. are possible, if stable under vacuum.
- Solid sample size:     (i) 5 mm (*l*) X 5 mm (*b*) X 2 mm (*h*)  
                                     (ii) 5 mm (*dia*) X 2 mm (*h*)
- Powder sample weight: 20 mg

## Types of analysis

- Survey Scan
- Narrow Scan
- Mapping
- Depth profiling
- Angle-resolved XPS

## Applications of XPS

- Survey the atoms, percentage of the elements, chemical state, and electronic state of the elements within a material
- Resolve mixtures of chemical states (e.g. Si, SiO, and SiO<sub>2</sub>)
- Profile the change in element atom% versus depth
- Map the uniformity of surface chemistry of a large area
- Surface contamination
- Chemical state changes or contamination due to processing or heating
- Check the purity of a chemical (e.g MoO<sub>3</sub>)
- Corrosion study
- Chemical depth profile down to 2,000 nm

### Details of XPS

Brand	Thermo Scientific
Model	K-Alpha
Sponsored Agency	RUSA 2.0 Physical Sciences - TRP (Ref. No. 63iii-1/BDU/RUSA2.0/TRP/PS Date:06.06.2022)



scan for more details.